Oxygen Calibration

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To perform the oxygen calibration, tap on \times until « O_2 Calibration» is shown on the screen. Choose with the O-button between a one-point calibration and two-point calibration and follow the instructions on the screen. When carrying out the two-point oxygen calibration, ensure you use 100% oxygen. For further information please refer to the user manual.

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Push and hold the ×-symbol for 5 seconds for the zero calibration.

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User Manual

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The user manual for the CITREX H4 can be found on our website at the link below. Please read the operating instructions carefully and use the product as recommended.

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Read the user manual carefully before using the CITREX H4. Online at www.imtanalytics.com/analyser_manuals

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
FlowLab Software FlowLab Software

View your flow, pressure and volume measurements comfortably on your computer screen while simultaneously saving respiratory parameters.

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Quick start manual - CITREX H4 Kurzanleitung - CITREX H4

Calibration

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Worldwide express 10-day service including shipping. More information on: www.imtanalytics.com/easycal

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To guarantee the accuracy of your measurements, your CITREX H4 must be recalibrated by IMT Analytics yearly.

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

Change Settings

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Tap the ×-symbol to open the settings menu. By tapping this symbol you can scroll through the settings. Use the O-symbol to change the

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14 RS-232

Numerische Werte

Tapping on the □-symbol will scroll through the available screens. By tapping on the symbol will scroll through the six configurable screens.

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Panels

Messkurven

Open and change the waveform screen by tapping the \triangle -symbol.

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Change parameter and units

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If the context button (≡) is pressed twice in quick succession, the "Edit Mode" is activated. A red icon on the screen indicates this. The parameter or the unit in the red frame can be changed with the \triangle -symbol or the \square -symbol. The O-symbol has the function to jump to the next parameter or unit. If the context button or the ×-button is pressed once, the "Edit Mode"

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Setup for Flow Measurement

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We recommend using a filter and the laminar flow tube for maximum accuracy. The filter protects the device from contamination and dust. Please make sure the filter is placed right before the laminar flow tube.

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
The measured gas must be free of oil, grease and dust.

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The device setup can affect the accuracy of the flow measurement. Tight bends, kinks or dents should be avoided. For more information please consult the user manual.

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Bad setup: Kink, tees, angle pieces at the device inlet

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Measurement Setup for Mechanical Ventilators

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The CITREX H4 is a powerful device for ventilator calibration and verification. We recommend you follow the measurement setup below.

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IMT Analytics AG . Gewerbestrasse 8 . 9470 Buchs . Switzerland T +41 81 750 67 10 . www.imtanalytics.com