

FlowAnalyser set VA

Biomedical test set

analyser
the art of measuring

Check and calibrate all types of gas flow and pressure-producing devices in a very short time. Create your individual report for each test.

The biomedical test set FlowAnalyser set VA is a certified, high precise measuring system delivered with all needed component parts.

Biomedical test set FlowAnalyser set VA to test and verify:

- Ventilators CPAP/Bilevel
- Ventilators ICU
- Ventilators Infant
- Ventilators High Frequency
- Anesthesia machines
- Anaesthetic vaporizer
- Blood pressure analysers
- Oxygen concentrators
- Vaccum pumps
- Spirometers
- Pipe gases

FlowAnalyser PF-300

The FlowAnalyser measures flow, pressure, temperature, humidity and O₂ concentrations bi-directionally. The device for measurements and calibration combines a simple, intuitive multilingual user interface with the highest precision.

FlowLab Software

FlowLab is the ideal software solution which has an impressive simple menu navigation.

Biomedical test set FlowAnalyser set VA includes:

FlowAnalyser PF-300/PF-301/PF-302, MultiGasAnalyser OR-703, SmartLung Adult, FlowLab Software, adapterset, bacterial filter, USB cable, power supply unit, carrying case, user manual

Order number: 300.116.004



IMT Analytics

FlowAnalyser set V

Biomedical test set

analyser
the art of measuring

Check and calibrate all types of gas flow and pressure-producing devices in a very short time. Create your individual report for each test.

The biomedical test set FlowAnalyser set V is a certified, high precise measuring system delivered with all needed component parts.

Biomedical test set FlowAnalyser set V to test and verify:

- Ventilators CPAP/Bilevel
- Ventilators ICU
- Ventilators Infant
- Ventilators High Frequency
- Blood pressure analysers
- Oxygen concentrators
- Vaccum pumps
- Spirometers
- Pipe gases

FlowAnalyser PF-300

The FlowAnalyser measures flow, pressure, temperature, humidity and O₂ concentrations bi-directionally. The device for measurements and calibration combines a simple, intuitive multilingual user interface with the highest precision.

FlowLab Software

FlowLab is the ideal software solution which has an impressive simple menu navigation.

Biomedical test set FlowAnalyser set V includes:

FlowAnalyser PF-300/PF-301/PF-302, SmartLung Adult, FlowLab Software, adapter-set, bacterial filter, USB cable, power supply unit, carrying case, user manual

Order number: 300.116.003



IMT Analytics

IMT Analytics AG . Gewerbestrasse 8 . 9470 Buchs . Switzerland
T +41 81 750 67 10 . www.imtanalytics.com